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
- Avoid eating, drinking alcohol, smoking, or exercising for at least 30 minutes before taking a measurement.
- Do not move or talk during measurement.

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Correct Posture

Sit upright with your back straight.

Remove tight fitting clothing from your upper arm, and any thick clothing such as a sweater. Do not place the cuff over thick clothes and do not roll up your sleeve so that it becomes tight around the arm.

Place your arm on a table so that the cuff will be at the same level as your heart.

The gap between the chair and the top of the table should be between 25cm to 30cm.

Common Causes of Inaccurate Results

Arm cuff is too loose.


Movement or talking during measurement.

Leaning forwards.

Note: If the arm cuff is at a lower position than your heart use cushions etc., to adjust the height of your arm.

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


Maintenance

To protect your unit from damage, please observe the following:

- Do not subject the main unit and cuff to extreme temperatures, humidity, moisture or direct sunlight.
- Do not fold the cuff or tubing tightly.
- Do not inflate the arm cuff over 299 mmHg.
- Do not disassemble the unit.
- Do not subject the unit to strong shocks or vibrations (for example, dropping the unit on the floor).
- Do not use volatile liquids to clean the main unit.
- Do not wash the arm cuff or immerse it in water.
- Do not use petrol, thinners or similar solvents to clean the arm cuff.
- Do not carry out repairs of any kind by yourself. If a defect occurs, consult your OMRON retail outlet or distributor as mentioned on the packaging.


- The unit should be cleaned with a soft, dry cloth.
- Use a soft, moistened cloth and soap to clean the arm cuff.

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8. Technical Data

Product Description	Digital Automatic Blood Pressure Monitor
Model	OMRON M2 Basic (HEM-7116-E/HEM-7116-E8)
Display	LCD Digital Display
Measurement Method	Oscillometric method
Measurement Range	Pressure: 0 mmHg to 299 mmHg Pulse: 40 to 180/min.
Memory	Last measurement
Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of display reading
Inflation	Automatic by electric pump
Deflation	Automatic pressure release valve
Pressure Detection	Capacitive pressure sensor
Power Source	4 "AAA" batteries 1.5V or AC/ DC adapter (optional, 6V $\approx$ 4W)
Battery Life	Capacity of new alkaline batteries is approx. 300 measurements
Operating Temperature/ Humidity	10°C to 40°C/ 30 to 85% RH
Storage Temperature/ Humidity/ Air Pressure	-20°C to 60°C/ 10 to 95% RH/ 700-1060 hPa
Console Weight	Approximately 245g without batteries
Cuff Weight	Approximately 120g
Outer Dimensions	Approximately 104 (w) mm $\times$ 64 (h) mm $\times$ 129 (l) mm
Cuff Dimensions	Approximately 146 mm $\times$ 446 mm (Medium cuff: arm circumference 22 to 32 cm)
Cuff material	Nylon
Package Content	Main unit, medium cuff, instruction manual, storage case, battery set, guarantee card, blood pressure pass

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

None-invasive blood pressure monitor

Display blood pressure(mmHg)

SYS = Systolic / DIA = Diastolic

Arm cuff

Display pulse rate
(BPM) ❤️
86 BPM


Important:

- If your systolic or diastolic pressure is outside the standard range, the heartbeat symbol will blink when the measurement result is displayed.

Recent research suggests that the following values can be used as a guide to high blood pressure for measurements taken at home.

Systolic Blood Pressure	Above 135 mmHg
Diastolic Blood Pressure	Above 85 mmHg

This criteria is for home blood pressure measurement. For professional office blood pressure measurement criteria, please refer to Chapter 9 "Some Useful Information about Blood Pressure".

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ECRI


Comparison

PROPOSED DEVICE INFORMATION		REQUIREMENTS	BIDDER'S RESPONSE		
Brand/Manufacturer Name	Bidder to specify		Nihon Kohden	Philips Healthcare	ZOLL Medical
Model Name	Bidder to specify		AED-3100	HeartStart FRx	AED Plus
Country of Origin	Bidder to specify		Japan	USA	UK
Year of 1st Production for Model Proposed	Bidder to specify		Not specify	2005	2002
FDA Clearance	Yes		JP FDA (Not sold in USA)	Yes	Yes
CE Mark (MDD)	Yes		Yes	Yes	Yes
Picture of Model Proposed	Bidder to attach a picture of model proposed				
CLINICAL APPLICATIONS/PROCEDURES		REQUIRED SPECIFICATIONS	BIDDER'S RESPONSE		
For external defibrillation to the heart of patients to restore normal rhythm and contractile function in patients who are experiencing ventricular fibrillation (VF), ventricular tachycardia (VT), or another shockable rhythm.		For general emergency care use	Yes	Yes	Yes
For monitoring physiologic parameters for patients who are suffering of cardiac arrests.		For general emergency care use	Yes	Yes	Yes
TECHNICAL SPECIFICATIONS/FEATURES		REQUIRED SPECIFICATIONS	BIDDER'S RESPONSE		
Type	Suitable for adult and pediatric patients		Yes	Yes	Yes
Configuration	Mounted on Crash Cart		Mounted on AED Box	Mounted on AED Box	Mounted on AED Box
Modes of Operation	3 modes		2 modes	2 modes	2 modes
a) External Defibrillation	Operator applies the paddles to patient's chest and discharges the defibrillator by simultaneously pressing two discharge buttons on paddle		Yes	Yes	Yes
Energy Selection, J	50J to 200J or 50J to 360J biphasic; For pediatric/neonatal: 2J to 20J		200 adult, 50 pediatric and infant	150 adult, 50 pediatric and infant	120, 150, 200 adult 50, 70, 85 pediatric

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

BiomedicalBenchmark™ THE TECHNOLOGY SUPPORT SYSTEM		
215 Categories; Updated Annually		
Device Description Universal Medical Device Nomenclature™	UMDNS™ Codes	Expected Life Years
Lights, Surgical (LED)	12-282	10
Lights, Ultraviolet, Germicidal	12-284	10
Linear Accelerators; Radiotherapy Units, Cobalt	16-972	12
Lithotripters, Extracorporeal	16-758	7
Lithotripters, Intracorporeal	16-230	8

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


Donate

Dispose/ Scrap

Note: Deletion of patient information

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Dropdown menu (highlighted with a red circle):

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